

AutoReceive System Requirements

Review the hardware and software requirements outlined below before you install AutoReceive. For optimal performance, maintain the preferred hardware requirements.

Microsoft® Windows Vista

Note: As of the AutoReceive 4.0 release, the Windows NT, Windows 98, and Windows Me operating systems are no longer supported.

Item	Requirement
Processor	<i>Preferred:</i> 1 GHz or higher <i>Minimum:</i> 800 MHz or higher
Memory	<i>Preferred:</i> 1 GB RAM <i>Minimum:</i> 512 MB RAM
Hard Disk	<i>Disk Preferred:</i> 2 GB with 400 MB of free space <i>Minimum:</i> 2 GB with 100 MB of free space
Monitor	800 x 600 through 1400 x 1050 resolution
Disk Drive	DVD ROM or CDROM drive
Printer	Laser printer connected to a workstation or office network
Internet Connection	56 KB minimum modem speed for ISP connection