

Medi-Gate Bridge Implementation

Installation Checklist

Use this form to ensure all steps required to install the Medi-Gate Bridge have been completed. If you need additional instruction, refer to the Medi-Gate Integrator™ Install Guide or the appropriate Quest Diagnostics bridge documentation.

	Install Step Description	Date Completed
<input type="checkbox"/>	Verify System Requirements on Upload Workstation <ul style="list-style-type: none"> • Operating System (Windows® 98, NT® 4.0 - SP 6a required, Windows ME, Windows 2000 – SP2 Recommended, XP) • Workstation Internet Connectivity 	
<input type="checkbox"/>	Verify Software Requirements on Upload Workstation <ul style="list-style-type: none"> • Microsoft Internet Explorer v. 5.01+ • Microsoft .NET Framework 	
<input type="checkbox"/>	Verify/Complete TOROL Administrative Setup (create bridge definitions)	
<input type="checkbox"/>	Install the Medi-Gate Integrator™ on Upload Workstation	
<input type="checkbox"/>	Configure the Medi-Gate Integrator™ <ul style="list-style-type: none"> • Enter OrgID, User, and Password 	
<input type="checkbox"/>	If PMS file creation method is “Print PMS report to Windows print driver/Medi-Gate printer”, then install the re-direction port monitor, the Medi-Gate/PostScript Print Drivers, and configure the printer port(s)	
<input type="checkbox"/>	If PMS file creation method is “Print PMS data to file”, then configure this method on the appropriate workstation plus the file transfer process (file must be accessible from Upload Workstation)	
<input type="checkbox"/>	Identify Appropriate Reports <ul style="list-style-type: none"> • Demographics Report Identified: _____ • Patient Insurance Report Identified: _____ OR <ul style="list-style-type: none"> • Name of Report Containing both: _____ 	
<input type="checkbox"/>	Train Customer how to Submit Reports	
<input type="checkbox"/>	Transmit Initial Report Information to Medi-Gate (for template creation)	
<input type="checkbox"/>	Complete Insurance Carrier Crosswalk	
<input type="checkbox"/>	Communicate Live Status (have customer submit first report for data population)	