

To determine the Operating System:

Right click the *My Computer* icon on the desktop
Go to *Properties> General Tab*
The operating system is noted under the heading: "System:"

To determine the Internet Explorer version:

Double click on your Internet Explorer icon (not necessary to be connected to Internet)
Go to *Help>About> Internet Explorer*
Note the Cipher Strength indicates the encryption level

To determine if Microsoft .NET Framework is installed:

Click on *Start, Settings, Control Panel*
Click on *Add/Remove Programs*
Scroll the list and check for the name

Workstation Requirements

Supported Operating Systems	Windows® 98 Windows NT® 4.0 (SP 6a required) Windows ME Windows 2000 – SP Recommended Windows XP
Communications	Internet Connectivity is required for each workstation that will upload patient data from the PMS
Browsers	Internet Explorer (IE) 5.01 or higher (Netscape not supported)
Other Software	Microsoft .NET Framework

PMS Requirements

- Refer to the PMS Qualification flowchart

TOROL Requirements

- Refer to the TOROL Qualification flowchart
- Refer to the TOROL Installation Scenario document