

Care360 Labs & Meds System Requirements

In order to access Care360 Lab & Meds, your workstation must meet the hardware and software requirements listed in the following table.

Item	Requirement
Processor	Minimum: Pentium II Preferred: Pentium III
Memory	Minimum: 128 MB Preferred: 512 MB <i>For Microsoft® Windows® XP Professional: 512 MB</i>
Hard disk	<i>Minimum: 100 MB of free hard disk space Preferred: 250 MB of free hard disk space</i>
Monitor	1024 x 768 resolution
Video Card	<i>Minimum: 256 colors Preferred: 16 bit color</i>
Printer	Laser printer connected to a workstation or office network
Label Writer	Dymo LabelWriter® XL for labels Zebra® Label printer model TLP3844-Z
Operating system	Microsoft® Windows XP Professional, Windows 2000, Windows Vista™, or Windows XP Pen Tablet PC Edition
Web browser	Internet Explorer 6.0, 7.0, & 8.0 with 128-bit encryption
Internet connection	Broadband; for example, Digital Subscriber Line (DSL) or cable modem
Browser Plug-ins	Adobe Reader 7.0 or newer Java Runtime Environment (JRE) 1.4.1 or newer Adobe Flash Player 9.0 or newer