

Technical Bulletin

Secondary Result Feed – Report Explanation

Issued: March 26, 2007

Summary:

The Secondary Result Feed Daily Report is sent via email between 8:30PM – 9:00 PM EST. The current distribution list consists of the Client Systems Managers. If you want changes made to who receives this in your BU, please email Sean Howard.

The daily report is based on the GMT time zone. When it is noon on the East coast, the servers see it as 4PM. This is important to consider when looking at dates. On the servers, the next day starts at 8PM EST.

The report only shows batches of lab results that have finished processing entirely.

Here is the process it goes through:

- A lab result has to come down via Lightning
- Copies over the network to LightningRouter
- Streams into the eGate (MedPlus) interface
- Imports into the Physician Portal.

Important to note that since these batches are larger than normal, the lag time is longer. For example, a single 9K lab result batch may take 3 hours just for the Lightning piece to finish. Queuing results in the afternoon – 2:00 PM (6:00 PM server time) – may not appear on that day's report.

Below is a sample of the daily report:

Care360 Physician Portal Secondary Feed Lab Activity
03/22/07 00:00 - 23:59:59 GMT

BU	Account	Count	Activity
AHL	22683752	18	Import successful
AHL	22683752	26	Newer version of message exists
DAL	18932	12	Import successful
DAL	48920	3	Import successful
HOU	3040300	65	Import successful
NOR	41722000	1	Import successful
QBA	236042	10	Import successful
QBA	268059	3686	Import successful
QBA	373083	121	Import successful
QBA	373083	1	Newer version of message exists
QBA	385463	1	Import successful
QDV	00301	8	Import successful

QDV 03908	1 Import successful
QDV 06092	1 Import successful
QDV 06117	2 Import successful
QDV 06118	2 Import successful
QDV 06119	3 Import successful
QDV 06120	3 Import successful
QDV 06540	1 Import successful
QDV 07142	2 Import successful
QDV 07159	3 Import successful
QDV 07164	5 Import successful
QDV 07230	6 Import successful
QDV 07245	1 Import successful
QDV 07266	3 Import successful
QDV 07268	5 Import successful
QDV 07277	1 Import successful
QDV 07291	2 Import successful
QDV 07295	2 Import successful
QDV 07350	6 Import successful
QDV 07366	10 Import successful
QDV 07370	1 Import successful
QDV 07384	6 Import successful
QPT 41447	55 Import successful
QPT 41448	28 Import successful
QPT 42114	6 Import successful
QPT 42115	10 Import successful
QPT 42183	18 Import successful